

made of biomedical stainless steel

developed tooth anatomy

*according to Dr. Washbon's study
in a wide selection of sizes*

graduated thickness

by variation in tooth size

"uniform stress" system shaping

*optimal elastic properties
over the entire area of the band*

laser etching

*strong retention on the tooth.
Ideal for REP, Herbst,
First Class appliances*

laser etched sizing

*it allows the sterilization
at high temperatures and
a durable good view of sizing*

free prewelding service

*Customized prewelding to any tube,
bracket or accessory of the section "G".
Ready delivery of standard prewelded bands
for the techniques most in use*



wide choice

of assortments in autoclavable trays

